



Recombinant *Oryza sativa* subsp. japonica Vacuolar iron transporter homolog 3 (Os08g0153300, LOC_Os08g05720)

Product Code	CSB-CF771178OFG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q84ZM7
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Sequence	MAMQMNSVHVHSTSPSPSPATSPPEEGKQEHGEVAAVHVVGVDDEAVMVV KDEEAFGGG GVDYSGRAQWLRAAVLGANDGLVSVASLMIGVGAVSESGRAMLVSGVAGLV AGAC SMAIG EFVSVYAQYDIEVAAARRRRRQRRRRCDGDGEEEGSGRLPSPFKAAAASALA FTVGALLP LLAGGFVRPWAPRVAAVCAATSAALAGFGALGAALGGASPARSAARVLLGGW AAMAACYG VLRLFANLY
Source	in vitro E.coli expression system
Target Names	Os08g0153300
Protein Names	Recommended name: Vacuolar iron transporter homolog 3 Alternative name(s): Protein NODULIN-LIKE 3
Expression Region	1-249
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.