



Recombinant Schizosaccharomyces pombe

Protein YIP5 (yip5)

Product Code	CSB-CF892549SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9UTD3
Product Type	Transmembrane Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Sequence	MAHEIEIDSEADLGRTTFEADLYKQRTPIPKPPAPRLQSRRASMQETPWTIKD SFNVET KDVVQRCIHTLIPTVNFFDVVDDRPDLYGPFWITTTVIQALFFSNSITEYARYAT GHGTS GYSIKKLISAASIIYGYTTIIAVLLWGILVWNKCNPKLLDCLCLYGYANIVWLPVSL ATP PFGLLSTLASHIVKYVLTGIGLLISIVFLTRNLYPICQQAGSNLCKLLLFGIIVFHCL LA LSLQLIFFS
Source	in vitro E.coli expression system
Target Names	yip5
Protein Names	Recommended name: Protein YIP5 Alternative name(s): YPT-interacting protein 5
Expression Region	1-249
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.