



Recombinant Human Membrane-spanning 4-domains subfamily A member 8B (MS4A8B)

Product Code	CSB-CF883629HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9BY19
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MNSM TSAVPVANSVLVVAPHNGYPVTPGIMSHVPLYPNSQPQVHLVPGNPPS LVS NVNGQ PVQKALKEGKTLGAIQIIIGLAHIGLSIMATVLVGEYLSISFYGGFPFWGGLWFII SGS LSVAAENQPYSYCLLSGSLGLNIVSAICSAVGVILFITDLSIPHPYAYPDYYPYA WGVNP GMAISGVLLVFCLEFGIACASSHFGCQLVCCQSSNVSVIYPNIYAANPVITPEP VTSPPSYSSEIQANK
Source	in vitro E.coli expression system
Target Names	MS4A8
Protein Names	Recommended name: Membrane-spanning 4-domains subfamily A member 8B Alternative name(s): Four-span transmembrane protein 4
Expression Region	1-250
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.