



Recombinant Clc-like protein 2 (clc-2)

Product Code	CSB-CF612085CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q11085
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MSQAVSYAILVLTIIAFLLTAAALCTPAWQVVYAREIRQWVQSGLWLSCQTRPN GMYSCT YTFSHDDFNTYFSDEVSGFRTPSFYPWQRTLFIHYLISQAFAMLSLISFCVSVS HKESKM PNILRSVFLVLAAVIAFGCLIAFAVYSYMVEYRFFHVSVSGIYEKHRGYSWYIAL TGAFV YLVAILSVVHVLLQARNSTTMSRQINSSLQSDFFEYQYHPNRSMESFEDRF AMRTL P PVPRQEKKTTVF
Source	in vitro E.coli expression system
Target Names	clc-2
Protein Names	Recommended name: Clc-like protein 2
Expression Region	1-252
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.