



Recombinant Bovine CD151 antigen (CD151)

Product Code	CSB-CF675839BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q3ZBH3
Product Type	Transmembrane Protein
Immunogen Species	Bos taurus (Bovine)
Sequence	MGEFGKSTTCGTVCLKYLLFTFNCCFWLAGLAVMAVGIWTLALKSDYISLLAS GTYLAT AYILVVAGIVVMVTGALGCCATFKERRNLLRLYFGLLLIIFLLEIIAGALAYIYYQQ LNA ELKENLKDTMTRRYHQPGHEGVTSAVDKLQQEFHCCGSNNSRDWQDSEWIH SGEAGGRVV PDSCCKTVVPGCGRRDHASNIYKVEGGCITKLETFIQEHLRIIGAVGLGIACVQV FGMLF TCCLYKSLKLEHY
Source	in vitro E.coli expression system
Target Names	CD151
Protein Names	Recommended name: CD151 antigen Alternative name(s): CD_antigen= CD151
Expression Region	1-253
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.