



Recombinant Mouse BRI3-binding protein (Bri3bp)

Product Code	CSB-CF805986MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8BXV2
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MGARASQEPRTVRAGLRVLLPVLLLALLLALVAPGAQGARGRGAADKNSH RRATSSFS QSVSSLFGEDNVRAAQKLLSRLTERFVQGVD MFLET LWKVWMELEVLGLDV SNLSQYFS PASVSNPTRLVLVGVVLLAYWFLSLTLGFTFSLHLVFGRRFFWLVRVILFSM SCVYIL HKYEGEPEHAVLPLCVVVAIYFMTGPMGYWRGSPGGLCSPSVEEKLEHLENQ VRLLNIRLNRVLENLDRSKDK
Source	in vitro E.coli expression system
Target Names	Bri3bp
Protein Names	Recommended name: BRI3-binding protein Short name= I3-binding protein
Expression Region	1-253
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.