



Recombinant Human Tetraspanin-11 (TSPAN11)

Product Code	CSB-CF025148HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A1L157
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MAHYKTEQDDWLIYLYKYLFFVFNFFFWGGAAVLAVGIWTLVEKSGYLSVLAS STFAAS AYILIFAGVLMVMTGFLGFGAILWERKGCLSTYFCLLLVIFLVELVAGVLAHVYY QRLSD ELKQHNLNRTLAENYGQPGATQITASVDRLQQDFKCCGSNSSADWQHSTYILL REAEGRQV PDSCKTVVVRCGQRAHPSNIYKVEGGCLTKLEQFLADHLLLMGAVGIGVACL QICGMVL TCCLHQRLQRHFY
Source	in vitro E.coli expression system
Target Names	TSPAN11
Protein Names	Recommended name: Tetraspanin-11 Short name= Tspan-11
Expression Region	1-253
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.