



Recombinant Chlorocebus aethiops CD151 antigen (CD151)

Product Code	CSB-CF004880DSU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P61170
Product Type	Transmembrane Protein
Immunogen Species	Chlorocebus aethiops (Green monkey) (Cercopithecus aethiops)
Sequence	MGEFNEKKTTCGTVCLKYLLFTYNCCFWLAGLAVMAVGIWTLALKSDYISLLAS GTYLAT AYILVVAGAVVMVTGVLGCCATFKERRNLLRLYFILLIIFLLEIIAGVLAYVYYQQ LNT ELKENLKDTMAKRYHQPGHEAVTSAVDQLQQEFHCCGSNNSQDWRDSEWIR LREARGRVV PDSCCKTVVAGCGQRDHAFNIYKVEGGFITKLETFIQEHLRVIGAVGTGIACVQ VFGMIF TCCLYRSLKLEHY
Source	in vitro E.coli expression system
Target Names	CD151
Protein Names	Recommended name: CD151 antigen Alternative name(s): Platelet-endothelial tetraspan antigen 3 Short name= PETA-3 CD_antigen= CD151
Expression Region	1-253
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.