



Recombinant Human Protein orai-2 (ORAI2)

Product Code	CSB-CF853504HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q96SN7
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MSAELNVPIDPSAPACPEPGHKGMDYRDWVRRSYLELVTSNHHSVQALSWR KLYLSRAKL KASSRTSALLSGFAMVAMVEVQLETQYQYPRPLLIAFSACTTVLVAVHLFALLIS TCILP NVEAVSNIHNLNSISESPHERMHPYIELAWGFSTVLGILLFLAEVLLCWIKFLP VDARR QPGPPPGPGSHTGWQAALVSTIIMVPVGLIFVVFTIHFYRSLVRHKTERHNREI EELHKL KVQLDGHHERSLQVL
Source	in vitro E.coli expression system
Target Names	ORAI2
Protein Names	Recommended name: Protein orai-2 Alternative name(s): CAP-binding protein complex-interacting protein 2 Transmembrane protein 142B
Expression Region	1-254
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.