



Recombinant Chicken Protein orai-2 (ORAI2)

Product Code	CSB-CF720030CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5ZLW2
Product Type	Transmembrane Protein
Immunogen Species	Gallus gallus (Chicken)
Sequence	MSSELNVPVDPSTPACCSEPGTKGMDYRDWVRRSYLELVTSNHHSVQALSW RKLYLSRAK LKASSRTSALLSGFAMVAMVEVQLEMQYKYPQMLLIAFSACTTVLVAVHLFALL ISTCIL PNVEAVSNIHNLNSISESPHERMHPYIELAWGFSTVLGILLFLAEVVLLCWIKFL PVDSI LKNDDTTTIQKPSGHAGWQAALVSTIIMVPVGLIFVVFTIHFYRSLVRHKTERHNR EIEEL HKLKVQLDGHDRGMQVV
Source	in vitro E.coli expression system
Target Names	ORAI2
Protein Names	Recommended name: Protein orai-2 Alternative name(s): Transmembrane protein 142B
Expression Region	1-257
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.