



Recombinant Mouse Uroplakin-1a (Upk1a)

Product Code	CSB-CF871951MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9D132
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MASAATEGEGKSPVVVGLLVVGNIIILLSGALFAETVWVTADQYRVYPLMGVS GKDDVF AGAWIAIFCGFSFFVVASFGVGAALCRRRYMILTYLLLMLIVYIFECASCITSYTH RDYM VSNPSLITKQMLTYYSADTDQGQELTRLWDRIMIEQECCGTSGPMDWVNYTS AFRAATPE VVFPWPPLCCRRTGNFIPINEDGCRVGHMDYLFTKGC FEHIGHAIDSYTWGIS WFGFAIL MWTLPVMLIAMYFYTTL
Source	in vitro E.coli expression system
Target Names	Upk1a
Protein Names	Recommended name: Uroplakin-1a Short name= UP1a Alternative name(s): Uroplakin Ia Short name= UPIa Short name= UPKa
Expression Region	1-257
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.