



Recombinant Bovine Melatonin receptor type 1A (MTNR1A)

Product Code	CSB-CF015189BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O02769
Product Type	Transmembrane Protein
Immunogen Species	Bos taurus (Bovine)
Sequence	YPLALASIVNDGWSLSSLHCQLSGFLMGLSVIGSVFNITGIAINRYCCICHSLRY NKLYS STNSLCYVFLIWMLTLVAIVPNLCVGTQLQYDPRIYSCTFTQSVSSAYTIAVVVFH FIVPM LVVIFCYLRIWALVLQVRWRVKPDNPKLKPQDFRNFMFVVFVLFVFAICWAPL NFIGLV VASEPASMARIPEWLFVASYMGYFNSCLNIIYGLLNQNFQRQEYRKIIIVSLCT TKMFF VDSSNHVAHRIKRKPS
Source	in vitro E.coli expression system
Target Names	MTNR1A
Protein Names	Recommended name: Melatonin receptor type 1A Short name= Mel-1A-R Short name= Mel1a receptor
Expression Region	1-257
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.