



Recombinant Human Uroplakin-1a (UPK1A)

Product Code	CSB-CF025654HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O00322
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MASAAAAEAEKGGSPVVVGLLVVGNIIILLSGLSLFAETIWVTADQYRVYPLMGV SGKDDV FAGAWIAIFCGFSFFMVASFGVGAALCRRRSMVLTYLVMLIVYIFECASCITSY THRDY MVSNP SLITKQMLTFYSADTDQGQELTRLWDRVMIEQECCGTSGPMDWVNFT SAFRAATP EVPFPWPPLCCRRTGNFIPLNEEGCRLGHMDYLFTKGCFEHIGHAIDSYTWGI SWFGFAI LMWTLPVMLIAMYFYTML
Source	in vitro E.coli expression system
Target Names	UPK1A
Protein Names	Recommended name: Uroplakin-1a Short name= UP1a Alternative name(s): Tetraspanin-21 Short name= Tspan-21 Uroplakin 1a Short name= UPIa Short name= UPKa
Expression Region	1-258
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.