



Recombinant Danio rerio Transmembrane protein 41A-A (tmem41aa)

Product Code	CSB-CF681374DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q502G2
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	VFLPPGPQLHKQSHEGETTDAKDGDEPSEMETASSRLKFPSDLDELKEMAEL LQFYKTEH TGYVLLLFC SAYLYKQAF AIPGSSFLN ILAGALFGTWFG LLLTCVLT TVGATLCF LLSQA FGKHHIVKLFDPKV AMLQKKVEENRSSLFFFLLFLRFFP MSPNWFLNMTSPILNI PVTLF FMAVFIGLMPYNFICVQTGSMLSQISSLDDLFSWSVVLKLLL TACVALLPGALIR KYSTR HLHLDGLETNGLSQNKKNR
Source	in vitro E.coli expression system
Target Names	tmem41aa
Protein Names	Recommended name: Transmembrane protein 41A-A
Expression Region	23-281
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.