



Recombinant Human Protein unc-50 homolog (UNC50)

Product Code	CSB-CF687496HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q53HI1
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MLPSTSVNSLVQGNGVLNSRDAARHTAGAKRYKYLRRLLFRFRQMDFEFAAW QMLYLFTSP QRVYRNFHYRKQTKDQWARDPAFLVLLSIWLCVSTIGFGFVLDMGFFETIKLL LWVLI DCVGVGLLIATLMWFISNKYLVKRQSRDYDVEWGYAFDVHLNAFYPLLVLHFI QLFFIN HVILTDTFIGYLVGNTLWLVAVGYYIYVTFGLGYSALPFLKNTVILLYPFAPLILLYG LSL ALGWNFTHTLCSFYKYRVK
Source	in vitro E.coli expression system
Target Names	UNC50
Protein Names	Recommended name: Protein unc-50 homolog Alternative name(s): Periodontal ligament-specific protein 22 Short name= PDLs22 Protein GMH1 homolog Short name= hGMH1 Uncoordinated-like protein
Expression Region	1-259
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.