



Recombinant Human Uroplakin-1b (UPK1B)

Product Code	CSB-CF025655HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O75841
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Transmembrane Proteins
Immunogen Species	Homo sapiens (Human)
Sequence	AKDNSTVRCFQGLLIFGNVIIGCCGIALTAECIFFVSDQHSLYPLLEATDNDDIYG AAWI GIFVGICLFCLSVLGIVGIMKSSRKILLAYFILMFIVYAFEVASCITAATQQDFFTP NLF LKQMLERYQNNSPPNDDQWKNNGVTKTWDRMLQDNCCGVNGPSDWQK YTSAFRTENND ADYPWPRQCCVMNNLKEPLNLEACKLGVPGFYHNQGCYELISGPMNRHAWG VAWFGFAIL CWTFWVLLGTMFYWSRIEY
Source	in vitro E.coli expression system
Target Names	UPK1B
Protein Names	Recommended name: Uroplakin-1b Short name= UP1bAlternative name(s): Tetraspanin-20 Short name= Tspan-20 Uroplakin 1b Short name= UPIb
Expression Region	1-260
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.