



Recombinant Arabidopsis thaliana Derlin-1 (DER1)

Product Code	CSB-CF819823DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8VZU9
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MSSPGEFYNSLPPITKAYGTLCCFFTTVATQLGLVAPVHIALIPELVLKQFQIWRLI TNLF FLGGFSINFGIRLLMIARYGVQLEKGPFFERRTADFLWMMIFGSFTLLVLSVIPFF WTPFL GVSLVFMLLYLWSREFPNANISLYGLVTLKAFYLPWAMLALDVIFGSPIMPDLL GIIAGH LYYFLTVLHPLATGKNYLKTPKWVNKIVARWRIGAPVASVRQAGGVGAAGPGA GGGVGGG GAYSSARAPPESSNTAFRGRSYRLTD
Source	in vitro E.coli expression system
Target Names	DER1
Protein Names	Recommended name: Derlin-1 Alternative name(s): AtDerlin1-1
Expression Region	1-266
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.