



Recombinant Human Probable palmitoyltransferase ZDHHC12 (ZDHHC12)

Product Code	CSB-CF822219HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q96GR4
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MAPWALLSPGVLVRTGHTVLTWGITLVFLHDTELQWEEQGELLPLTFLLLV LGSLLL YLAVSLMDPGYVNVQPPQEQEELKEEQTAMVPPAIPLRRCRYCLVLQPLRARH CRECRRCV RRYDHHCPWMENCVGERNHPLFVVYLALQLVLLWGLYLAWGLRFFQPWG QWLRSSGLL FATFLLLSLFLVASLLLVSHLYLVASNTTTFEFISSHRIAYLRQRPSNPFDRGL TRNLA HFFCGWPSGSWETLWAEIEEEEGSSPAV
Source	in vitro E.coli expression system
Target Names	ZDHHC12
Protein Names	Recommended name: Probable palmitoyltransferase ZDHHC12 EC= 2.3.1.- Alternative name(s): Zinc finger DHHC domain-containing protein 12 Short name= DHHC-12 Zinc finger protein 400
Expression Region	1-267
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.