



Recombinant Human Membrane-spanning 4-domains subfamily A member 12 (MS4A12)

Product Code	CSB-CF885756HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9NXJ0
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MMSSKPTSHAEVNETIPNPYPPSSFMAPGFQQPLGSINLENQAQGAQRAQPY GITSPGIF ASSQPGQGNIQMINPSVGTAVMNFKEEAKALGVIQIMVGLMHIGFIVLCLISFS FREVL GFASTAVIGGYPFWGGLSFIISGSLSVSASKELSRCLVKGSLGMNIVSSILAFIG VILLL VDMCINGVAGQDYWAVLSGKGISATLMIFSLLLEFFVACATAHFANQANTTTNM SVLVIPN MYESNPVTPASSSAPPRCNNYSANAPK
Source	in vitro E.coli expression system
Target Names	MS4A12
Protein Names	Recommended name: Membrane-spanning 4-domains subfamily A member 12
Expression Region	1-267
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.