



Recombinant Mouse Gap junction gamma-3 protein (Gjc3)

Product Code	CSB-CF842087MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q921C1
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MLLLELPIKCRMCGRFLRQLLAQESQHSTPVGRFLLPMLMGFRLLILVSSGPGV FGNDEN EFICHLGQPGCKTICYDVFRPLSPLRFWAFQVILMAVPSAIYVAFTLYHVIGYWE VPGKE NKEQETQISKGDHSGKDVSGAKSLKLLWAYVAHLGVRLALEGAALGVQYNLYG FKMSSTFI CREDPCIGSTTCFQSHPSSEKTIFLNIMFGISGACFLFIFLELALLGLGRFWRIYKH KLSF LKKLPTSESSVRSKDTTDELSVVEAKEPF
Source	in vitro E.coli expression system
Target Names	Gjc3
Protein Names	Recommended name: Gap junction gamma-3 protein Alternative name(s): Connexin-29 Short name= Cx29 Gap junction epsilon-1 protein
Expression Region	1-269
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.