



Recombinant Danio rerio Transmembrane protein 192 (tmem192)

Product Code	CSB-CF023784DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6NYE7
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MVMESGRTSAYGNNTNVDFQSTDDDDPLIDGPLVSQDALESAIKREFQTLSTH WAVGVLS FLHVYVVVVCIVITVFCWISDEHTSECTAALNGIDSKTVVLLGKVALWVVVFIYD RVVQR RHSVRRKGYLDFYRMTKGIKSLPLLIHSAGNAAILTVIAPSSLLDANVKNLSVY LLAI ICLELLLTIIICLLRYTVHVVKFNRKKPHPDVTEEERSYGCSSDAHTETGFRDGSS LEDVV EKQADLIDYLKQHNSVLSRRILTLTGQQIRA
Source	in vitro E.coli expression system
Target Names	tmem192
Protein Names	Recommended name: Transmembrane protein 192
Expression Region	1-271
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.