



Recombinant Human Probable lipid phosphate phosphatase PPAPDC3 (PPAPDC3)

Product Code	CSB-CF818757HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NBV4
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MPASQSRARARDRNNVLNRAEFLSLNQPPKGGPEPRSSGRKASGPSAQPPP AGDGARERR QSQQLPEEDCMQLNPSFKGIAFNSSLAIDICMSKRLGVCAGRAASWASARSM VKLIGITG HGIPWIGGTILCLVKSSTLAGQEVLNLLLALLLDIMTVAGVQKLIKRRGPYETS PSLLD YLTMDIYAFPAGHASRAAMVSKFFLSHLVLAVPLRVLLVLWALCVGLSRVMIGR HHVTDV LSGFVIGYLQFRLVELVWMPSSSTCQMLISAW
Source	in vitro E.coli expression system
Target Names	PLPP7
Protein Names	Recommended name: Probable lipid phosphate phosphatase PPAPDC3 EC=3.1.3.- Alternative name(s): Phosphatidic acid phosphatase type 2 domain-containing protein 3
Expression Region	1-271
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.