



Recombinant *Saccharomyces cerevisiae* Vacuolar membrane protein YPL162C (YPL162C)

Product Code	CSB-CF604896SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q12042
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Sequence	MYVSN GK DTCQLLGPVSLFVQTLMGMTAVIVLLVKRNYEHPRRKMIWWSYDIG KQIIIGSL GIHFLNLGISILKKRRRSLFAITAKGNDDDEDQCDWYFLNLLDTTVGIPILWLCLY IIEK VLKSLHFQNIESGNYFPSKTVGSHPRKPLFSAFVKQLLIFIVGLGVMKFCVFLIL NYLED LAYWFADLILGWSDSWPNFQVFLVMFVFPILLNCFQYFCVDNVIRLHSESLTIT NAENFE TNTFLNDEIPDLSEVSNEVPNKDNNISSYGSII
Source	in vitro E.coli expression system
Target Names	YPL162C
Protein Names	Recommended name: Vacuolar membrane protein YPL162C
Expression Region	1-273
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.