



Recombinant Human E3 ubiquitin-protein ligase MARCH5 (MARCH5)

Product Code	CSB-CF889137HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NX47
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MPDQALQQMLDRSCWVCFATDEDDRTAEWVRPCRCRGSTKWWHQACLQR WVDEKQQRGNST ARVACPQCNAEYLIVFPKLGPVVYVLDLADRLISKACPF AA AGIMVGSIIWTAV TYGAVT VMQVVGHKGLDVMERADPLFLLIGLPTIPVMLILGKMIRWEDYVLRRLWRKYS NKLQILN SIFPGIGCPVPRIPAEANPLADHVSATRILCGALVFPTIATIVGKLMFSSVNSNLQ RTIL GGIAFVAIKGAFKVYFKQQQYLRQAHRKILNYPEQEEA
Source	in vitro E.coli expression system
Target Names	MARCH5
Protein Names	Recommended name: E3 ubiquitin-protein ligase MARCH5 EC= 6.3.2.- Alternative name(s): Membrane-associated RING finger protein 5 Membrane-associated RING-CH protein V Short name= MARCH-V Mitochondrial ubiquitin ligase Shor
Expression Region	1-278
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.