



Recombinant Squalus acanthias Proteolipid protein DM beta

Product Code	CSB-CF336244FKW
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P36964
Product Type	Transmembrane Protein
Immunogen Species	Squalus acanthias (Spiny dogfish)
Sequence	MEENMEEGQTQKGCFECCIKCLGGIPYASLIATILLYAGVALFCGCGHEALSGT VTILQN NFEVVRGAGDTIDVFTMIDIFKYVIYGVAIAFFVYGILLMVEGFFTTGAIKDLYGD FKIT TCGRCVSGWFIMLTIFMLAWLGVTAFSTLPVFMFYFNIWTLQCQNVTIMESTDLC FDLRQF GIVPIHEQKTVCTLNENFSKLCQSNLDNMTFHLFIVALAGAGAAVIAMVHYLMV LSANWA YVKDACRMQKYEDIKSKEEQELHDIHSTRSKERLNAYT
Source	in vitro E.coli expression system
Protein Names	Recommended name: Proteolipid protein DM beta
Expression Region	1-278
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.