



Recombinant Human Gap junction gamma-3 protein (GJC3)

Product Code	CSB-CF854118HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NFK1
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MCGRFLRLLAEESRRSTPVGRLLLLPVLLGFRLVLLAASGPGVYGDEQSEFVC HTQQPGC KAACFDAFHPLSPLRFWVFQVILVAVPSALYMGFTLYHVIWHWELSGKGKEEE TLIQGRE GNTDVPGAGSLRLLWAYVAQLGARLVLEGAALGLQYHLYGFQMPSSFACRRE PCLGSITC NLSRPSEKTIFLKTMTFGVSGFCLLFTFLELVLLGLGRWWRTWKHKSSSSSKYFL TSESTRR HKKATDSLPVVETKEQFQEAVPGRSLAQEKQRPVGPRDA
Source	in vitro E.coli expression system
Target Names	GJC3
Protein Names	Recommended name: Gap junction gamma-3 protein Alternative name(s): Connexin-30.2 Short name= Cx30.2 Connexin-31.3 Short name= Cx31.3 Gap junction epsilon-1 protein
Expression Region	1-279
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.