



Recombinant Mouse Leukocyte antigen CD37 (Cd37)

Product Code	CSB-CF713974MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q61470
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MSAQESCLSLIKYFLFVFNLFVGLGGLIFCFGTWILIDKTSFVSVFVGLSFVPLQT WSKV LAVSGVLTMALALLGCVGALKELRCLLGLYFGMLLLLFATQITLGILISTQRVRLE RRVQ ELVLRTIQSYRTNPDETAEEESWDYAQFQLRCCGWQSPRDWNKAQMLKANE SEEPFVPCS CYNSTATNDSTVFDKLFSSQLSRLGPRAKLRQTADICALPAKAHIYREGCAQSL QKWLHN NIISIVGICLGVGLLELGFMTLSIFLCRNLDHVYDRLARYR
Source	in vitro E.coli expression system
Target Names	Cd37
Protein Names	Recommended name: Leukocyte antigen CD37 Alternative name(s): CD_antigen= CD37
Expression Region	1-281
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.