



Recombinant Arabidopsis thaliana Secretory carrier-associated membrane protein 1 (SCAMP1)

Product Code	CSB-CF882895DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9SKT3
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MSRYQSHSFDDGEINPFANPTSVPAATSKLSPLPEPYDRGATMDIPLDSGKD LKAKEKE LREKEAELKRREQEIKRKEDAIAQAGIVIEEKNWPPFFPLIHHDISNEIPIHLQRIQ YVA FTSMLGLVVCLLWNIVAVTTAWIKGEGPTIWFLAIIYFISGVPGAYVMWYRPLYR AMRTD SALKFGWFFFTYLFHIAFCVFAAVAPPIIFKGKSLTGILPAIDVLSGNILVGIFYFIG FG FFCLESLVSIWVIQQVYMYFRGSGKAAEMKQEATRRAMMAAL
Source	in vitro E.coli expression system
Target Names	SCAMP1
Protein Names	Recommended name: Secretory carrier-associated membrane protein 1 Short name= AtSC1 Short name= Secretory carrier membrane protein 1
Expression Region	1-282
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.