



Recombinant Danio rerio Presqualene diphosphate phosphatase (ppapdc2)

Product Code	CSB-CF713149DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5TZ07
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MPSPKARSGSGRSGSVPCPGGNGRYEFISLNRTPPSPVPPQLLQRQGS DPTT ARLRASES PVRRRGSGSSNSSTGGGGQQLPEEDCMRLNPSFFGIALSSLLAIDLWLSKRL GVCACEDS SWGSVRPLMKLIEVSGHGIPWLAGAAYCLYKSDSPAGQEVMLNLLMALVLDV VLVGVLKA VVRRRRPAHNRMDFATFSVDSYSFPSGHATRAAMCARFLNHLVLAAPLRV LVLLWATI VGFSRVLLGRHNVTDVAFGFFMGYWQYNLVEMLWLSPVMLQSAIGQLH
Source	in vitro E.coli expression system
Target Names	p1pp6
Protein Names	Recommended name: Presqualene diphosphate phosphatase EC= 3.1.3.- Alternative name(s): Phosphatidic acid phosphatase type 2 domain-containing protein 2
Expression Region	1-288
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.