



Recombinant Arabidopsis thaliana Secretory carrier-associated membrane protein 3 (SCAMP3)

Product Code	CSB-CF873464DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9M5P2
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MANRYDPNPFEEEEVNPFANARGVPPASNSRLSPLPPEPVGFDYGRTVDIPL DRAGTQD LKKKEKELQAKEAELKRREQDLKRKEDAAARAGIVIEVKNWPPFFPLIHHDIAN EIPVHL QRLQYVTFATYLGLVLCLFWNIIAVTTAWIKGEGVTIWLLALIYFIAGVPGGYVL WYRPL YRAFRTDSALSFGWFFLFYMLHIAFCVF AAVAPPVVFKGKSLAGILPAIDVLSG QAIVGI FYFIGFAFFCLESVSVSIWVIQQVYMYFRGSGKQDQMRREAARGALRAAV
Source	in vitro E.coli expression system
Target Names	SCAMP3
Protein Names	Recommended name: Secretory carrier-associated membrane protein 3 Short name= AtSC3 Short name= Secretory carrier membrane protein 3
Expression Region	1-289
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.