



Recombinant Mouse Protein orai-3 (Orai3)

Product Code	CSB-CF750827MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6P8G8
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MKGGEGDTGEQAPLNPEVDSPAGSATYREFVHRGYLDLMGASQHSRLRALSW RRLYLSRAK LKASSRTSALLSGFAMVAMVEVQLENDHEYPPGLLVAFSACTTVLVAVHLFAL MVSTCLL PHIEAVSNIHNLNSVHQSPHQRLHRYVELAWGFSTALGTFLFLAEVVLVGWVK FVPIGAP MGKPAPVVPMSQVPPVTVSLSLASNLTPSSASITTSQQPSKACPPRQVCDSA HGPGWQAA MASTAIMVPVGLVFMALHFYRSLVAHKTDNRHKQELEELSRLQGELQAV
Source	in vitro E.coli expression system
Target Names	Orai3
Protein Names	Recommended name: Protein orai-3 Alternative name(s): Transmembrane protein 142C
Expression Region	1-290
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.