



Recombinant Bacillus subtilis Uncharacterized transporter ydbO (ydbO)

Product Code	CSB-CF308732BRJ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P96610
Product Type	Transmembrane Protein
Immunogen Species	Bacillus subtilis (strain 168)
Sequence	MERTENLKKGEK GALLNIFAYVILAVV KLVIGILYHSEALRADGLNNGTDIVASVA VLIG LRISQRPADSDHPYGHYRAETIS SLVASFIMMAVGIEVLIGGGKAIAGGTTETPN LIAAW TALGSAVFMYGIIYLYNKRLAASIKSSALMAAAKDSRSDAFVSAGAFIGVFSSQL KLPWVD PVTAFIIGIIICKTAWDIFK DASHSLTDGFHLKDLEPYKQTVGRIENVHRLKDVKA RYLG STVHIEMVITVDPKLTVEEGHGV ADEVEDKIKHEHDVTHVHVHVEPDDIK
Source	in vitro E.coli expression system
Target Names	ydbO
Protein Names	Recommended name: Uncharacterized transporter ydbO
Expression Region	1-290
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.