



Recombinant Prolipoprotein diacylglyceryl transferase (Igt)

Product Code	CSB-CF352734EGX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P60957
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Sequence	MTSSYLHFPEFDPVIFSIGPVALHWYGLMYLVGFIFAMWLATRRANRPGSGWT KNEVENL LYAGFLGVFLGGRIGYVLFYNFPQFMADPLYLFRVWDGGMSFHGGLIGVIVVM IIFARRT KRSFFQVSDFIAPLIPFGLGAGRLGNFINGELWGRVDPNFPFAMLFPGSRTEDI LLLQTN PQWQSIFDTYGVLPRHPSQLYELLLEGVVLFIILNLYIRKPRPMGAVSGLFLIGY GAFRI IVEFFRQPDAQFTGAWVQYISMGQILSIPMIVAGVIMMVWAYRRSPQQHVS
Source	in vitro E.coli expression system
Target Names	Igt
Protein Names	Recommended name: Prolipoprotein diacylglyceryl transferase EC= 2.4.99.-
Expression Region	1-291
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.