



Recombinant Mouse B-lymphocyte antigen CD20 (Ms4a1)

Product Code	CSB-CF015007MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P19437
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MSGPFPAEPTKGPLAMQPAPKVNLKRTSSLVGPTQSFFMRESKALGAVQIMN GLFHITLG GLLMIPTGVFAPICLSVWYPLWGGIMYIISGSLIAAAAETSRKSLVKAKVIMSS LSLFA AISGIILSIMDILNMTLSHFLKMRRLELIQTSKPYVDIYDCEPSNSSEKNSPSTQY CNSI QSVFLGILSAMLISAFFQKLVTAGIVENEWKRMCTRSKSNVLLSAGEKNEQTI KMKEEI IELSGVSSQPKNEEEEIEIIPVQEEEEEEAEINFPAPPQEQESLPVENEIAP
Source	in vitro E.coli expression system
Target Names	Ms4a1
Protein Names	Recommended name: B-lymphocyte antigen CD20 Alternative name(s): B-cell differentiation antigen Ly-44 Lymphocyte antigen 44 Membrane-spanning 4-domains subfamily A member 1 CD_antigen= CD20
Expression Region	1-291
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	10xHis-SUMO-tag
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.