



Recombinant Mouse Protein YIF1A (Yif1a)

Product Code	CSB-CF842039MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q91XB7
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MAYHSAYGVHGSKHRTRAAPDPPPLFDDTSGGYSSQLGGYPAPGADVAFSV NLLGDPVA NMAMAYGTSIASQGKDIVHKELHRFVSVNKLKYFFAVDTAYVAKKLGLLVFPYT HQNWKM QYSHDVPLPPRKDLNAPDLYIPTMAFITYVLLAGMALGIQQRFSPEVLGLCAST ALVWVF MEVLALLLGLYLATVRSELSTFHLLAYSGYKYVGMILSVLTGLLFGSDGYVAL AWTSSA LMYFIVRSLRTAASGPDSMGGPAPRQRLQLYLTLGAAAFQPLIIYWLTFFHLVR
Source	in vitro E.coli expression system
Target Names	Yif1a
Protein Names	Recommended name: Protein YIF1A Alternative name(s): YIP1-interacting factor homolog A
Expression Region	1-293
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.