



Recombinant Human Gap junction delta-3 protein (GJD3)

Product Code	CSB-CF847609HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8N144
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MGEWAFGLSLLDAVQLQSPLVGRLWLVMILIFRILVLATVGGAVFEDEQEEFV CNTLQPG CRQTCYDRAFPVSHYRFWLFHILLSSAPPVLFVVYSMHRAGKEAGGAEAAAQ CAPGLPEA QCAPCALRARRARRCYLLSVALRLLAELTFLGGQALLYGFRVAPHFACAGPPC PHTVDCF VSRPTEKTVFVLFYFAVGLLSALLSVAELGHLLWKGRPRAGERDNRCNRAHEE AQKLLPP PPPPPPPALPSRRPGPEPCAPPAYAHPAPASLRECGSGRGKASPATGRRDL AI
Source	in vitro E.coli expression system
Target Names	GJD3
Protein Names	Recommended name: Gap junction delta-3 protein Alternative name(s): Connexin-31.9 Short name= Cx31.9 Gap junction alpha-11 protein Gap junction chi-1 protein
Expression Region	1-294
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.