



Recombinant *Ashbya gossypii* Protein ARV1 (ARV1)

Product Code	CSB-CF755097DOT
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q750Q7
Product Type	Transmembrane Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Sequence	MICINCCCHVESLYVAYPGDHIRLTDCWQCGEVVDRYVEFDNVLLFIDLLLLKP GAYRHL VYNSLEVDMRRYPEHGTRESGFRGYVQDWLLWFRKYDRLNRIWLLLITFEVYL MWMTEEK KYNYYREPATQYGWHLTGDVLSLQNPLLQYLFFAVYVIGDLALLHKL TRESL LNWSQW GQDLRYAKDIVSYTILLSYGAKIFPILMLIWPYDSL VSTSIKWIANFYIVESLRIVT SK SYSYIILLFAGIFIVRSVVLKPV LALIVSRGNMLQAQRYLWGEYELFKFRFLTKRD IFL
Source	in vitro E.coli expression system
Target Names	ARV1
Protein Names	Recommended name: Protein ARV1
Expression Region	1-299
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.