



Recombinant Protease HtpX homolog (htpX)

Product Code	CSB-CF835918FMR
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q93D93
Product Type	Transmembrane Protein
Immunogen Species	Streptococcus mutans serotype c (strain ATCC 700610 / UA159)
Sequence	MLFNQIASNKRKTIVLLIVFFMLLAAIGA AVGYLWLD SLVGGMAIALIIGFIYAFSM IFQ STNIVMAMNNAKEITVDEMPDYHIVEDMAMVAQIPMPRVYVIDDPSLNAFATG SSPQNA AVAATTGLLAIMNREELEGVIGHEVSHIRNYDIRISTIAVALASAVTLISSIGGRM MWFG GGRSRRSNDNEGGGYLQIILLIFSLAILAPLAASLVQLAISRQREYLADASSVE LTRN PQGMIRALQKLEQSQPMASVDNASAALYINDPKKKSGLKHLFYTHPPIADRIR RLEQV
Source	in vitro E.coli expression system
Target Names	htpX
Protein Names	Recommended name: Protease HtpX homolog EC= 3.4.24.-
Expression Region	1-299
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.