



Recombinant *Methylophilum infernorum* Protease HtpX homolog (htpX)

Product Code	CSB-CF457109MSV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B3DX74
Product Type	Transmembrane Protein
Immunogen Species	<i>Methylophilum infernorum</i> (isolate V4) (<i>Methylophilum infernorum</i> (strain V4))
Sequence	MFKRIFLLTLTNIIVIFLLTLFISLLHLDRWLNAYGIDYQTLFLFSMVVGFTGSFIS LAI SKWMAKMAYNIHVIQEPSNEAERWLIVETVAELAKRANIRMPEVGIYESPEVNA FATGPSR SNALVAVSTGILSQMNKKQIAGVLAHEITHINNGDMVTMTLLQGVVNTFVVFLS RIIGFF VDRLFSTRNEERESIGIGFYLGMFISEIVLGLLASIIVAWYSRMREFRADAGGAHL AGKEA MISALKKQIMEGESAFIDERSPALNAFKINGRPGGILALLATHPPLDERIKALE RIPD
Source	in vitro E.coli expression system
Target Names	htpX
Protein Names	Recommended name: Protease HtpX homolog EC= 3.4.24.-
Expression Region	1-300
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.