



Recombinant Cholesterol 25-hydroxylase-like protein (F35C8.5)

Product Code	CSB-CF631321CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q20027
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MLDLYPVQNLTVDQLEYEKNTFLQPAWDWIKNGNEHILSSPLFPPFYALSIDY TWVAVF TFIDVFLCNVPPFFKDAKIQKDRKVTWDLIKKSLKLGWNQLLWIYPMALVQLIW VPDEL PILAPTVFEMLSQLAIFFLAFDFTYFWFHYINHVKWLYRWCHSVHHMYSSPFA ASQHL HPFELFFVGTFITTIPWIFPTHCLTYWIWFFIAQSVSYEYVHIGYDFPFALHRIFWF YSGA PAHDMHHLRPLTCFQPWFNYLDRLMGYHITYADLKKMTEAKFKKFGLYSAED EKGLIKIN
Source	in vitro E.coli expression system
Target Names	F35C8.5
Protein Names	Recommended name: Cholesterol 25-hydroxylase-like protein EC= 1.14.99.-
Expression Region	1-300
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.