



Recombinant Zea mays CASP-like protein 15

Product Code	CSB-CF476051ZAX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B6TWJ1
Product Type	Transmembrane Protein
Immunogen Species	Zea mays (Maize)
Sequence	MALQAQQQPTSPTRDRVGSGEWLADTEKLPGAAASPEDVVVASTHHAAAA ARYVPPRAT SHTAEPNPGRDGGGGWYSWNGARRARHDPPAPRRQQPAKTHPPAPPLPAA PPPPPPPPA ASPAPAPRAPPHPAQVRSADRVVPAILSRKRRAAVMQRAALLARAAAAGXXLA ALTVLAA DTRRGWARDSYSNYAQFRYSEAVNVVGFLYSVFQFVALAELMRRNTHLIPHP KRGLFDFT MDQVLAYLLISSSSSATARASDLTENWGSDSFPNMANGSIAISFVAFVVFAICS LISAYN LFRRDM
Source	in vitro E.coli expression system
Protein Names	Recommended name: CASP-like protein 15 Alternative name(s): UL36 tegument protein
Expression Region	1-306
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.