



Recombinant Human Protein YIPF1 (YIPF1)

Product Code	CSB-CF896523HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9Y548
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MAAVDDLQFEFEGNAATSLTANPDATTVNIEDPGETPKHQPGSPRGSGREED DELLGNDD SDKTELLAGQKKSSPFWTFEYYQTFDQVDTYQVFDRIKGSLLPIPGKNFVRLYI RSNPDL YGPFWICATLVFAIAISGNLSNFLIHLGEKTYHYVPEFRKVSIAATIIYAYAWLVPL ALW GFLMWRNSKVMNIVSYSFLEIVCVYGYSLFYIPTAILWIIPQKAVRWILVMIALGI SGS LLAMTFWPAVREDNRRVALATIVTIVLLHMLLSVGCLAYFFDAPEMDHLPTTTA TPNQTV AAAKSS
Source	in vitro E.coli expression system
Target Names	YIPF1
Protein Names	Recommended name: Protein YIPF1 Alternative name(s): YIP1 family member 1
Expression Region	1-306
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.