



Recombinant Human Tumor suppressor candidate 3 (TUSC3)

Product Code	CSB-CF025352HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q13454
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	QKKKENLLAEKVEQLMEWSSRRSIFRMNGDKFRKFIKAPPRNYSMIVMFTALQ PQRQCSV CRQANEEYQILANSWRYSSAFCKLFFSMVDYDEGTDVFFQQLNMNSAPTFM HFPPKGRPK RADTFDLQRIGFAAEQLAKWIADRTDVHIRVFRPPNYSGTIALALLVSLVGGLLY LRRNN LEFIYNKTGWAMVSLCIVFAMTSGQMWNHIRGPPYAHKNPHNGQVSYIHGSS QAQFVAES HIILVLNAAITMGMVLLNEAATSKGDVVGKRRICLVGLGLVFFFSFLLSIFRSKY HGYP YSDLDFE
Source	in vitro E.coli expression system
Target Names	TUSC3
Protein Names	Recommended name: Tumor suppressor candidate 3 Alternative name(s): Magnesium uptake/transporter TUSC3 Protein N33
Expression Region	42-348
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.