



Recombinant Danio rerio Protein FAM26E (fam26e)

Product Code	CSB-CF008198DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A4FUN9
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MDSFKTVLKFFTDQKTTIGYSVMAILTIGSERIFSMVSFQCPCTKGQNFYPYGICF LLGPA VLLLVGFFVSTRFWRLYTGCCLNPLKLCPRGNFVGCLKGLIKVLYGACVAPV MWLTVL LNGTFYECVAVSGLDEVAVIQIFCADKGSVCQDELHRVPCGKSTLPPEQNTTELL YMLRAQS QILGWSVIITAVVIALIGTCYKNCVSQVSYLQLTFWKIYMEKEREKFDADFANDYA TKLAD RNLKSFFDNKLPEEFPPFNHKAWEEISAIYNFRKSEQHYSTLQRYVERSDRDY SPDEHPM VEMDHGIEMV
Source	in vitro E.coli expression system
Target Names	calhm5.1
Protein Names	Recommended name: Protein FAM26E
Expression Region	1-310
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.