



Recombinant Synechocystis sp. Uncharacterized transporter sll1263 (sll1263)

Product Code	CSB-CF305007SSQ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P74068
Product Type	Transmembrane Protein
Immunogen Species	Synechocystis sp. (strain PCC 6803 / Kazusa)
Sequence	MTARLARPYAVLSIGAALATMGLKLGAYAITGSVGLLSDALESTVNLASAIIVAF WALSLA ATPADSEHPFGHskaeyfSSGLEGAFIFVAALGIGYSAVERLLSPRPLDQNALG IALAIA ATALNGTVAWILWRAGKRLNSIALRADSQHLMTDVWTSVGVVVAVALIFVTGW EWLDPLI ALGVGFNVLWTGTHLLRETISSLMDQSLPPAQLQAITSCFLPYEDQGVRFHLL QTRQAGS QSFISFHVLVPGHWTVQRGHDLCEAIETAIAERITGSRVTTHLEPLEDPKSWQH PDEFPP SAPLNRDKPN
Source	in vitro E.coli expression system
Target Names	sll1263
Protein Names	Recommended name: Uncharacterized transporter sll1263
Expression Region	1-310
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.