



Recombinant Schizosaccharomyces pombe Palmitoyltransferase PFA5 (pfa-5)

Product Code	CSB-CF863184SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9C0W9
Product Type	Transmembrane Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Sequence	MSTNLEEKQRRFSKLQKTLGVVFPAAILLSTGYTVVWFIALICVDSNIKIRNGYR NLGGG IVLIFFFFITSGLAYFSYFRVLFSSPSFCGNTLYTYYGFDNPIFLCGPNGAPRMCG TCKC WLPDRSHHSRVSMRCIRKFDHYCSFVGKDVCFNSQKFFYQFLFYGFSAACMV LISTAIMI SRTYHYRSLPGTWIFVLVFSAFGVFLGVMLVRHTGYLLLNINSHEAKNWKTRI YSFSVF FPEHMDSRVLVQSDPGDLPWDRGYSENWRAVMGDHWYNWILPLRRSPGDG EHFLYSPSFV SKMQSKALSMNS
Source	in vitro E.coli expression system
Target Names	pfa-5
Protein Names	Recommended name: Palmitoyltransferase PFA5 EC= 2.3.1.- Alternative name(s): Protein fatty acyltransferase 5
Expression Region	1-312
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.