



Recombinant Human Protein FAM26F (FAM26F)

Product Code	CSB-CF722248HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5R3K3
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MEKFAVLDLHVKHHSALGYGLVTLTLAGGERIFSAVAFQCPCSAAWNLPYGL VFLLVPA LALFLLGYVLSARTWRLTGCCSSARASCGSALRGSLVCTQISAAAALAPLTW VAVALLG GAFYECAATGSAAFAQRLCLGRNRSCAAELPLVPCNQAKASDVQDLLKDLKA QSQVLGWI LIAVVIILLIFTSVTRCLSPVSFLQLKFWKIYLEQEQQILKSKATEHATELAKENIK CF FEGSHPKEYNTPSMKEWQQISSLYTFNPKGQYYSM LHKYVNRKEKTHSIRST EGD TVIPV LGFVDSSGINSTPEL
Source	in vitro E.coli expression system
Target Names	CALHM6
Protein Names	Recommended name: Protein FAM26F
Expression Region	1-315
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.