



Recombinant Human Protein YIPF2 (YIPF2)

Product Code	CSB-CF887141HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9BWQ6
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MASADELTFHEFEEATNLLADTPDAATTSRSDQLTPQGHVAVAVGSGGSYGA EDEVEEES DKAALLQEQQQQQPGFWTFSYYQSFFDVDTSQVLDRIKGSELLPRPGHNFVR HHLRNRPD LYGPFWICATLAFVLAVTGNLTLVLAQRRDPSIHYSPOFHKVTVAGISICYAWL VPLAL WGFLRWRKGVQERMGPYTFLETVCIYGYSLFVFIPMVVLWLPVWLPWLQWLFQ ALALGLSA AGLVFTLWPVVREDTRLVATVLLSVVLLHALLAMGCKLYFFQSLPPENVAPPP QITSLP SNIALSPTLPQSLAPS
Source	in vitro E.coli expression system
Target Names	YIPF2
Protein Names	Recommended name: Protein YIPF2 Alternative name(s): YIP1 family member 2
Expression Region	1-316
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.