



# Recombinant Human Zinc transporter 2 (SLC30A2)

|                          |                                                                                                                                                                                                                                                                                                                                                                        |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-CF874791HU                                                                                                                                                                                                                                                                                                                                                         |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.                                                                                                                                                                                                                                                                                                      |
| <b>Uniprot No.</b>       | Q9BRI3                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Product Type</b>      | Transmembrane Protein                                                                                                                                                                                                                                                                                                                                                  |
| <b>Immunogen Species</b> | Homo sapiens (Human)                                                                                                                                                                                                                                                                                                                                                   |
| <b>Sequence</b>          | MEAKEKQHLLDARPAIRSYTGSLWQEGAGWIPLPRPGLDLQAIELAAQSNHHC<br>HAQKGPD<br>SHCDPKKGKAQRQLYVASAICLLFMIGEVEILGALVSVLSIWVVTGVLVYLAVE<br>RLISG<br>DYEIDGGTMLITSGCAVAVNIIMGLTLHQSGHGHSHGTTNQQEENPSVRAAFIH<br>VIGDFM<br>QSMGVLVAAYILYFKPEYKYVDPICTFVFSILVLGTTLTILRDVILVLMEGTPKGV<br>DFTA<br>VRDLLSVEGVEALHSLHIWALTVAQPVLSVHIAIAQNTDAQAVLKTASSRLQG<br>KFHFHT VTIQIEDYSEDMKDCQACQGPSD |
| <b>Source</b>            | in vitro E.coli expression system                                                                                                                                                                                                                                                                                                                                      |
| <b>Target Names</b>      | SLC30A2                                                                                                                                                                                                                                                                                                                                                                |
| <b>Protein Names</b>     | Recommended name: Zinc transporter 2 Short name= ZnT-2 Alternative name(s): Solute carrier family 30 member 2                                                                                                                                                                                                                                                          |
| <b>Expression Region</b> | 1-323                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                    |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                          |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                    |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                     |