



Recombinant Human Sideroflexin-3 (SFXN3)

Product Code	CSB-CF880119HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9BWM7
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MESKMGELPLDINIQEPRWDQSTFLGRARHFFTVTDPRNLLLSGAQLASRNI VQNYRAG VVTPGITEDQLWRAKYVYDSAFHPDTGEKVVLIGRMSAQVPMNMTITGCMLTF YRKTPTV VFWQWVNQSFNAIVNYSNRSGDTPITVRQLGTAYVSATTGAVATALGLKSLTK HLPPLVG RFVPFAAVAAANCINIPLMRQRELQVGIPVADEAGQRLGYSVTAAKQGIFQVVI SRICMA IPAMAIPPLIMDTLEKKDFLKRPPWLGAQLQVGLVGFCLVFATPLCCALFPQKS SIHISN LEPELRAQIHEQNPSVEVVYYNKGL
Source	in vitro E.coli expression system
Target Names	SFXN3
Protein Names	Recommended name: Sideroflexin-3
Expression Region	1-325
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.